## GARANTIE

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## WARRANTY

This appliance is guaranteed for 12 months from the date of purchase. During the period of the guarantee, the distributor will take responsibility for defects occur, the appliance will be repaired or exchanged as the distributor deems necessary. Damage arising from improper use or normal wear and defects, which do not materially affect the performance of the appliance, is not covered under this guarantee. The guarantee is not valid if unauthorized repairs have been attempted or if replacement parts used are not approved by the supplier. The guarantee is valid only if it is stamped and signed by the retailer at the time of the purchase. For all claims made during the period of the guarantee, send the appliance and the warranty card to your retailer. All other claims are excluded from the guarantee unless specifically prohibited by law.

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**IMPORTANT SAFETY INSTRUCTIONS**

- It is recommended to use the appliance on a electrical installation using a differential with nominal current loss which doesn't exceeds 30mA.

- If the supply cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid a hazard.

- Do not reach for an appliance that has fallen into water. Unplug it immediately.

- Do not use while bathing or showering.

- Do not place or store appliance where it can fall or be pulled into a tub or sink.

**- Do not place in or drop into water or other liquid.**

**- You can fix the parts only if they cannot drop into water.**

- Do not attempt to plug or unplug unit while feet are in the water.  
- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts or attachments.

- Close supervision is necessary when this appliance is used by, on, or near children, invalids or disabled persons.

- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the supplier; specifically any attachments not provided with the unit.

- Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to appropriated service center for examination and repair.

- This appliance is designed for household use only.

- The supply cord cannot be replaced. If the cord is damaged the appliance should be scrapped.

**SAVE THESE INSTRUCTIONS**

**Specifications**

Charging unit: AC 230V, AC (~) 50/60 Hz. 2,5 W.

Rechargeable power unit: Battery: Ni-MH 1/2AA 300mAh 2.4V, max. operation time approx. 20/25 min. when fully charged.

Notice: when the rechargeable cells are running empty the plaque remover automatically turns off.

**Connecting**

The charging unit can be connected to any voltage outlet -even to shaver sockets- between 120- 240 volts AC (-).

Place the rechargeable power unit, facing the front onto the charging base.

The red indicator will light up when the toothbrush is inserted on the plugged in charging base.

The toothbrush may warm up somewhat during recharging. This has no significance on the unit's operation. Over-charging is impossible.

**To maintain the toothbrush's full performance, it should always be stored on the charging base after use.**

**Preserving the rechargeable batteries**

To maximize the capacity and life of the rechargeable batteries, they should be fully charged for 16 hours before the first use and then fully discharged by using Then recharge to full capacity again

This discharging/charging procedure should be done every six months.

**How to use**

- Before using your Plaque Remover for the first time, charge the power unit at least 16 hours

- After the power unit has been fully charged, attach the attachable brush head to the power unit.

- Before doing so, you may fix one of the brush head onto the end of the brush head.

- The brush head must face the same direction as the switch.

- Before brushing your teeth, rinse your mouth thoroughly and wet the brush head.

- You may use regular toothpaste. To avoid splashing, switch the power unit on after guiding the brush head to your teeth. Then move the brush head along the chewing surfaces of your teeth.

- Guide it along the outsides and then the insides of your teeth. Be sure to get into the spaces between your teeth as well as along the gum line

- You do not need to press hard To thoroughly remove plaque, brush for at least 2 minutes.

- After cleaning your teeth, rinse your mouth thoroughly again. During the first few days of using, your gums may bleed slightly, particularly if they normally tend to do so. After a few days, the bleeding will generally stop. If your gums continue to bleed after 2 weeks, you should see your dentist.

- Clinical tests have shown that the best cleaning effect is achieved by gently moving the brush from tooth to tooth, without exerting pressure

- A change in the shape of the brush is normal after use. This is influenced by your individual cleaning technique, but does not affect the efficiency of the brush.

**WARNING**

Do not use your device in a damp atmosphere (sauna) or close to an electrical appliance connected to the mains and the ground, or even piping connected to the ground.

**Cleaning**

To maintain the function of the brush head, rinse thoroughly for several seconds with the power unit switched on under running water after use. Clean the brush head and the power unit by simply holding them under running water. Wipe the power unit dry. After cleaning place the power unit back on the charging base. Additionally, we recommend soaking the brush head in a denture cleaner (commercially available) every 2 to 4 weeks to keep it hygienically clean.

**Protecting the environment**

-Remove the built-in rechargeable battery when eventually discarding your appliance.

-You can also take the appliance to a Service Centre, The staff will be happy to help you.

-Before removing, the battery should be completely empty of charge, Let the appliance run until the motor stops,

-Only unscrew the bottom for battery removal. The toothbrush may not be watertight anymore after the bottom has been unscrewed,

-This product contains NiMH cells. In the interest of protecting the environment, please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at a Service Centre or at appropriate collection sites provided in your country.

	Recommendations: When disposing of your product at end of life cycle, please comply with national regulations. Please dispose the device in allocation provided for this purpose in order to ensure safe and environmentally-friendly destruction.
---	--

**Norms :**

**CE : -(2004/108/CE/EMC)  
-(2006/95/CE/LVD)**

**RoHS (2002/95/EC) :** The product is made with friendly environmentally materials.

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